

# Tabuľka veľkostí IMPLEO™

Kód produktu	Rozmery základne cm	Projekcia cm	Váha g
<b>IMP-MR: Impleo stredná projekcia</b> 			
IMP-MR 150	9.4	3.4	165
IMP-MR 180	10.0	3.5	190
IMP-MR 210	10.5	3.6	220
IMP-MR 240	11.0	3.9	255
IMP-MR 270	11.2	4.1	278
IMP-MR 300	11.6	4.2	300
IMP-MR 330	12.0	4.3	340
IMP-MR 360	12.3	4.4	365
IMP-MR 390	12.7	4.4	390
IMP-MR 420	13.1	4.6	430
IMP-MR 460	13.6	4.7	480
<b>IMP-HR: Impleo vysoká projekcia</b> 			
IMP-HR 150	8.4	4.0	150
IMP-HR 180	8.9	4.2	180
IMP-HR 210	9.5	4.2	210
IMP-HR 240	10.0	4.4	240
IMP-HR 270	10.2	4.8	270
IMP-HR 300	10.7	4.9	300
IMP-HR 330	11.2	4.9	330
IMP-HR 360	11.6	5.0	360
IMP-HR 390	11.8	5.1	375
IMP-HR 420	12.4	5.1	425
IMP-HR 460	12.7	5.4	460
IMP-HR 520	12.9	5.6	505
IMP-HR 560	13.3	5.7	550
<b>IMP-EHR: Impleo extra vysoká projekcia</b> 			
IMP-EHR 150	8.2	4.3	160
IMP-EHR 180	8.4	4.6	180
IMP-EHR 210	9.0	4.8	210
IMP-EHR 240	9.4	4.8	235
IMP-EHR 270	9.8	4.9	260
IMP-EHR 300	10.4	5.0	290
IMP-EHR 330	10.9	5.2	330
IMP-EHR 360	11.2	5.4	360
IMP-EHR 390	11.5	5.5	400
IMP-EHR 420	11.9	5.6	435
IMP-EHR 460	12.2	5.8	460
IMP-EHR 520	12.6	6.1	520
IMP-EHR 560	13.0	6.4	575
IMP-EHR 620	13.3	6.6	620

Rozmery implantátov sa môžu jemne líšiť.

Pre viac informácií o produktoch navštívte webstránku [www.nagor.sk](http://www.nagor.sk).

© Copyright 2013 Nagor Ltd., 129 Deerdykes View, Westfield Industrial Estate, Cumbernauld, Glasgow, G68 9HN, United Kingdom